How physicians\patients can remotely upload images to Corewell Health East

- 1. Verify the patient's email in Epic. If they do not have an email in Epic, direct the patient to sign up for *my*BeaumontChart. This will create an email in Epic. Once this has been completed, you can email the patient a link to upload images to AMBRA.
- 2. Email this link to the physician\patient: <u>beaumont.ambrahealth.com/share/beaumont_fhir</u>
- 3. The physician\patient can "Create an Account" or "Login" and follow the steps below.

Beaumont	
	ImageShare
	Share medical images with Beaumont
	1 Enter your email address
	2 Select your image CD or file
	3 Upload and share your images
	Enter your email address
	Continue
	Indications for Use Terms of Use Privacy Policy

4. If they have an AMBRA account, it will ask them to log in. if they do not, it will have them create an account.

1 Enter your email address
2 Select your image CD or file
3 Upload and share your images
Enter password
Usernastmail@beaumonr.org
* Password
Forgot your password?
Continue

Share medical images with Beaumont

ImageShare

5. This is what creating an AMBRA account looks like:

ImageShare

Share medical images with Beaumont

1 Enter your email address

- 2 Select your image CD or file
- 3 Upload and share your images

Enter your name and a password

* Email
CreateNewEmail@Gmail.com
* First name
* Last name
* Password
8 characters or longer, containing 2 numbers and 2 letters
* Confirm password
* I accept the Terms of Use
* I accept the Privacy Policy
* I accept the Indications for Use
Continue

6. The next step is locating the images on CD or File location.

Beaumont

1							-				
I	2	Y	٦.	~	2	0	~	h	0	r	0
ı				ч	u	C	$\mathbf{\cdot}$		ч		
					0						

Share medical images with Beaumont

1	Enter	your	email	addr	ess
2	Selec	t you	r imag	e CD	or file

3 Upload and share your images

Need help? Watch a short video to learn how to find the right file on your computer:

×

To begin the upload process, select

Choose File.

Choose your images

Select files for DICOM wrapping (Ctrl+click to select multiple)



7. Navigate to the location of the images.

Look <u>I</u> n:	DVD RW Drive (D:)
	Л

- 8. Select the studies to upload.
 - a. Single or multiple studies can be selected to upload.
 - Any studies already in AMBRA will show RED and will not be allowed to be uploaded again. If it needs to be uploaded again, contact: <u>CHEReleaseOfImages@corewellhealth.org</u>

~

ImageShare

Shai	ihare medical images with Beaumont						
1 Enter your email address 2 Select your image CD or file 3 Upload and share your images Choose your images to upload and share							
	Description	Images	Patient Name	DOB	Modality	Study Date	
	Duplicate study found						
	CHEST	1	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
	kNEEE	2	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
v	cHEST	1	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
	BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	
	BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	
	BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	
ls this fo	or a consult?						
selec							
Messag	je:						
]*lu	inderstand this study will I	be shared with Beaumon	t, I have authority to share this st	tudy.			

Upload Selected Studies Preview Selected Studies Cancel

9. Click "Upload Select Studies."

ImageShare Share medical images with Beaumont 1 Enter your email address 2 Select your image CD or file 3 Upload and share your images

Choose your images to upload and share

Description	Images	Patient Name	DOB	Modality	Study Date	
Duplicate study found						
CHEST	1	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
KNEEE	2	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
cHEST	1	Pie ^ Cherry	10-04-1985	CR	03-03-2023	
BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	
BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	
BREAST	1	Pie ^ Cherry	10-04-1985	MG	03-03-2023	

Is this for a consult?

select

Message:

I understand this study will be shared with Beaumont, I have outhority to share this study.

Upbled Selected Studies
Preview Selected Studies
Cancel

10. When the study is uploaded, it will appear in the "Activities" queue in AMBRA for Radiology Support approval. Once the documentation from the user has been completed, it will appear in the chart.

Beaumont			
	ImageShare		
	 Your images were uploaded and shared with Beaumont 		
	Upload another study	View your images	Upload a report
	Choose another CD or image file to upload	View images you have uploaded and shared	Optionally upload a medical report if one was included with your medical images
	Upload another study	View images	Choose File No file chosen
		Indications for Use Terms of Use Privacy	Policy 03-01-23 04:26:19 UTC +AMBRHEALTHSOLUTIONS0/5\$+7_V3.23.1.0R

11. For questions, please contact <u>CHEReleaseOfImages@corewellhealth.org</u> or your local imaging support team.